

# OptiBond™ Selection Guide

## LEGEND

✓	Recommended
✓	Secondary Recommendation
✗	Not Recommended
-	Not Applicable
•	Yes
•	No



INDICATIONS				
Direct	✓	✓	✓	✓
Indirect	✗ <sup>1</sup>	✓	✓	✓
Anterior	✓	✓	✓	✓
Posterior	✓	✓	✓	✓
SUBSTRATES				
Dentin/Enamel/Uncut Enamel	✓	✓	✓	✓
All Ceramics (Including Lithium Disilicate Ceramic)	✗ <sup>1</sup>	✓ <sup>2</sup>	✓	✓
Metal Based (Including Zirconia Based & Alumina Based)	✗ <sup>1</sup>	✓ <sup>3</sup>	✓	✓ <sup>3</sup>
Inlay/Onlays	✗ <sup>1</sup>	✓	✓	✓
Core Build Up	✓	✓	✓	✓
Cementation of Posts	✗ <sup>1</sup>	✓	✓	✓
Composite/Metal Repair	✓	✓	✓	✓
Porcelain Repair	✓	✓	✓	✓
Amalgam Sealing	✓	✓	✓	✓
AUXILIARY PRODUCTS NEEDED				
Gel Etchant	•	•	✓ <sup>4</sup>	✓ <sup>4</sup>
Silane	-	•	• <sup>5</sup>	• <sup>5</sup>
Metal Primer	-	• <sup>5</sup>	• <sup>5</sup>	• <sup>5</sup>
Total Dark Cure Mode Capability for Cementation	-	•	• <sup>6</sup>	• <sup>6</sup>
BONDING COMPATIBILITY				
Compatible with all resin cements & core build up	-	• <sup>3</sup>	• <sup>7</sup>	• <sup>3</sup>

## FOOTNOTES

1. N/A Due To High Film Thickness
2. Use With Silane
3. Use Only With NX3/Maxcem Elite Cement
4. Optional
5. Not required
6. With NX3/Maxcem Elite/Maxcem Elite Chroma
7. Compatible with any resin cements.