

Quick Facts



AlgiNot™ Impression

Dimensionally stable, up to 30 days, easy to read, multiple pours.

Alginate Impression

Already exhibiting shrinkage within 24 hours.

AlgiNot FS™
Alginate Alternative

No more hand mixing. No lengthy clean up. Now that's a change worth making.

When does AlgiNot add value compared to traditional Alginates?

- You have patients returning with a broken temporary, but now you can easily repour your AlgiNot impression.
- Your AlgiNot study model can be shown during esthetic case discussions with your patients.
- Offer additional services with AlgiNot such as night guards, mouthguards, sportsguards, and bleaching trays.
- Use AlgiNot when offering sleep disorder treatments or denture preparation services.
- Replace your patient's broken sportsguards/nightguards quickly.*

What clinical applications is AlgiNot best used for?

- » Temporary crowns
- » Orthodontic models and study models
- » Removable prosthetics
- » Surgical impressions and guides
- » Sleep disorder treatments
- » Nightguards, mouthguards, sportsguards, bleaching trays



AlgiNot vs. Alginates

AlgiNot Alginate Alternative	Alginates
Available in a Cartridge or Volume Mixer. No dust, no hand mixing, no clean-up	Is mixed by hand. Clean and Lube for easy bowl cleanup
Can be disinfected easily (Compatible with CaviCide™)	Should not be disinfected as it can cause distortion
Good accuracy because of its high dimensional stability up to 30 days	Results in less accurate casts due to low dimensional stability (shrinkage)
Can be poured multiple times *	Can only be poured once
Can be poured immediately or readily stored for convenient use	If not poured immediately, must be kept moist for storage
Quick! Need versatility? Pour and repour when you want and need.*	Easy case – same day pour. DONE.

AlgiNot™ rebounds from deformation far better than traditional alginate “Jeltrate® Alginate” and better than alginate alternative “Position™ Penta™”

Visit kerrdental.com/alginot for more information.

*AlgiNot's dimensional stability allows you to pour and even repour your model up to 30 days. (Internal data on file)